

## Application Note AN R522

# Forensic Analysis of Documents Using Raman Spectroscopy

Documents have been investigated to determine the feasibility of utilizing Raman and SERS Raman spectroscopy for the identification and characterization of inks on paper. Fluorescence reduction methods have been employed to facilitate the analysis by reducing the nascent fluorescence from paper and ink. Further, ink crossings were investigated to demonstrate that ink applied after creation of a document could be differentiated from the originally applied ink.

### Introduction

Raman spectroscopy has been shown to be important in the forensic analysis of trace evidence, controlled substances, and documents<sup>1</sup>. Document analysis can be broken down into a few categories of investigation, including authentication of documents of value, determination of document alteration, and association of evidence to a suspect or place of origin. Because Raman spectroscopy, like infrared spectroscopy accesses fundamental modes of vibration, it is highly specific for chemical identification. Raman lines of condensed phase samples are typically narrower than that of IR absorptions (or IR bands), further aiding in the discrimination of molecular species. It is worth noting that Raman analysis is non-destructive and allows in-situ analysis even when applied to investigations in aqueous solu-


tions. Infrared microscopy measurements could be limited at longer wavelengths, due to the use of photoconductive mercury-cadmium-telluride (MCT) detectors that have a long wavelength cutoff of about  $600\text{ cm}^{-1}$ . Modern Raman microscopes routinely allow the detection and identification of inorganic pigments that frequently have spectral signatures down to  $\sim 100\text{ cm}^{-1}$ .

There are a few obstacles to the successful Raman analysis of documents. First, paper frequently exhibits fluorescence upon visible excitation, where the fluorescence signal dwarfs the Raman signal. Secondly, many synthetic dyes found in textiles, inks, and paints have inks and pigments that also exhibit fluorescence. Third, many inks are dark in color and subsequently undergo transient heating upon impingement of the incident laser light. This can yield damage or alteration of the document in question. Fluorescence reduction or rejection methods must therefore be frequently employed to collect useful Raman spectra. The first method employed for this investigation is the use of a concave baseline correction method that eliminates the perturbing fluorescence spectrum that is superimposed on the desired Raman spectrum. The second method used was Surface Enhanced Raman Scattering (SERS). The SERS effect is characterized by an increase in the Raman intensity by many


orders of magnitude for species adsorbed on rough metal (usually silver) surfaces compared to that obtained from the same number of molecules in solution or the gas phase<sup>3,4</sup>. At the same time, proximity to the surface provides a non-radiative pathway for relaxation from the excited state, which successfully quenches fluorescence. Recent work in (SERS) has demonstrated that organic colorants present in inks, paints, and textile fibers can be easily identified from microscopic samples by treatment of the sample with silver nanoparticles prior to analysis<sup>5</sup>. This method has proven effective for the detection of specific dyes from samples as small as a one-millimeter section of a single silk fiber of fifty-micrometer diameter and on dyes used in textiles (on archaeological samples) even when severely degraded by burial<sup>6</sup>. Emerging methods for completely non-destructive SERS analysis of documents and textile fibers will enable the identification of inks and dyes while preserving the integrity of evidence such as textile fibers or documents.

### Application Examples


One of the areas of interest is the detection and characterization of ink applied after the creation of a document. This would be important in determining whether a document has been altered and if so which inks were applied first. Confocal Raman microscopy can be applied such that a depth profile can be performed on the ink in question to determine which occurred first. The depth resolution for confocal Raman microscopy is up to 1 micron. If ink is applied first and allowed to dry, subsequent application of ink would sit on top of the earlier deposited ink. If an ink applied later is composed of a different pigment, Raman would readily differentiate between the inks. Figure 1 shows an ink crossing for two different inks applied within a few seconds of each other. Figure 2 shows the raw Raman spectra collected from the ball point pen (top) and the z-grip pen (bottom). The Raman spectra were collected using a SENTERRA Raman microscope (Bruker Optics, Inc.) at 532 nm with a 1s integration in the confocal mode of operation at a resolution of  $3\text{ cm}^{-1}$ . The SENTERRA was mounted on a z-stage to allow access to large samples as shown in Figure 3. Those same Raman spectra are shown in Figure 4 after removal of the fluorescence by utilization of the concave baseline correction method. Prior to fluorescence removal the overlaying Raman spectra are barely observable. The corresponding Raman images are shown in Figure 5, where the  $1648\text{ cm}^{-1}$  was integrated for the z-grip image and the  $1146\text{ cm}^{-1}$  band for the ball point pen image. The ball point pen was readily identified as methyl violet and the z-grip ink as methyl violet with an additive. It is clearly evident that the z-grip pen was applied over the ball point pen.




Visible images of two crossing inks on a document. The close-up image on the right shows the mapping grid used in acquiring the Raman images.




Raw Raman spectra of z-grip pen (top, offset) and ball point pen (bottom) collected at 532 nm.




SENTERRA Raman microscope mounted on a z-stage.

**Figure 4**

Concave baseline corrected Raman spectra from the Ball point pen (top, offset) and the z-grip pen (bottom).

**Figure 6**


Depth profile of z-grip pen as a function of depth.

**Figure 5**

Visible image of pen ink crossing (top left), Raman chemical images of paper ( $1088\text{ cm}^{-1}$  integration) (upper right), z-grip pen ( $1647\text{ cm}^{-1}$  integration) (lower left), and image of the ball point pen ( $1146\text{ cm}^{-1}$  integration) (lower right).

We can extend the investigation further by bringing the document through the field of focus to depth resolve each of the inks at the junction point, as shown for the z-grip pen in Figure 6. The depth profile was conducted by stepping in 1 micron steps from the top surface focus position into the sample for a total of 100 microns. As expected the spectrum at the interface was composed of a combination of the spectra representative of both inks (data not shown). This example demonstrates the facility of confocal Raman microscopy for document analysis, so long as the fluorescence can be effectively removed.


Another area of forensic interest is the analysis of documents generated from inkjet printers. Figure 7 shows the visible images from such a document. The resultant Raman spectra also exhibited significant fluorescence, where concave baseline correction was also employed. The Raman images, shown in Figure 8, demonstrate the ability of Raman spectroscopy to be used for the characterization of documents via pigment identification.

**Figure 7**

Visible images of a document formed from ink jet deposition of the ink. The black color clearly contains more than one pigment.


Another document was examined with SERS to determine if the fluorescence could be effectively quenched as well as achieve a signal enhancement. This document had ink containing a phthalo blue and a methyl violet ink. The 488 nm spectrum shows some peaks of both over a fluorescent background. The quality is low due to the paper fluorescence, and to the methyl violet fluorescence. Going to 785 nm does not improve things much because, as the fluorescence of methyl violet decreases, phthalo blue becomes slightly fluorescent. In this case, a droplet of silver colloid suspension, prepared using a previously described method where microwave assisted reduction is employed, was deposited onto the document using a piezo device prior to analysis<sup>7,8</sup>. The silver colloid droplet size was about 20 microns in diameter. 488 nm excitation was used for the analysis of this document. Figure 9 shows the Raman spectrum without correction (top) and the SERS result (bottom). It is clear that the fluorescence has been effectively quenched and a very high quality SERS Raman obtained. This striking example shows the promise of using SERS Raman microscopy in the detection of very small amounts of dyes and pigments even on fluorescing substrates such as paper or even skin.

Figure 8



Visible image of ink jet document (top left), Raman chemical images of black color (1000 - 1400  $\text{cm}^{-1}$  region integration), red color (200-600  $\text{cm}^{-1}$  region integration) (lower left), and image of the blue color (750  $\text{cm}^{-1}$  peak integration) (lower right).

Figure 9



Raman spectrum of the ball point ink on paper with 488 nm excitation (top) and the corresponding SERS spectrum after silver colloid deposition (bottom).

## Conclusions

In summary, Raman microanalysis can become an important tool in the forensic analysis of documents. Previously, the interference from fluorescence precluded the widespread application of Raman to this field of study. With the advent of current state-of-the-art methods for removing fluorescence and the emergence of SERS Raman as a viable and practical method for quenching fluorescence and enhancing the Raman signal, it is expected that interest and activity in this field will increase greatly.

<sup>1</sup> Bartick, E.G., Applications of Vibrational Spectroscopy in Criminal Forensic Analysis, in Handbook of Vibrational Spectroscopy, J.M. Chalmers and P.R. Griffiths (Eds), John Wiley and Sons Ltd, Chichester, 2002.

<sup>2</sup> M. Pirzer and J. Sawatzki, US Patent 2006/0212275.

<sup>3</sup> Lombardi, J.R.; Birke, R.L.; Lu, T.; Xu, J., Charge-transfer theory of surface enhanced Raman spectroscopy; Herzberg-Teller contributions; J. Chem. Phys. 1986, 84, 4174.

<sup>4</sup> Rodger C, Dent G, Watkinson J, Smith WE., Surface Enhanced Resonance Raman Scattering and Near-Infrared Fourier Transform Raman Scattering as in situ Probes of Ink Jet Dyes Printed on Paper; Applied Spectroscopy 54, 11, 2000 1567-1576.

<sup>5</sup> Review Papers, 14th International Forensic Science Symposium, Interpol-Lyon, 19-22 October 2004, N. N. Daeid (Editor).

<sup>6</sup> Leona M, Stenger J, Ferloni E. Application of surface-enhanced Raman scattering techniques to the ultra-sensitive analysis of natural dyes in works of art. ArtRaman 2005 3rd International Conference on the Application of Raman Spectroscopy in Art and Archaeology, Paris (F), 31 August 2005 - 3 September 2005.

<sup>7</sup> Leona, M, Microanalysis of organic pigments and glazes in polychrome works of art by surface-enhanced resonance Raman scattering. Proc Natl Acad Sci USA 106:14757-14762, 2009.

<sup>8</sup> Leona, M, Tague Jr, TJ, US Patent Pending.

**Thomas J. Tague Jr., Ph.D. and Peng Wang, Ph.D.**

Bruker Optics, Inc.,  
19 Fortune Drive, Billerica, MA 01821

**Ph.D. Marco Leona,**

Metropolitan Museum of Art,  
1000 Fifth Avenue, New York, NY 10028

### ● Bruker Optik GmbH

Ettlingen · Deutschland  
Phone +49 (7243) 504-2000  
Fax +49 (7243) 504-2050  
info@brukeroptics.de

[www.bruker.de/optik](http://www.bruker.de/optik)

### Bruker Optics Inc.

Billerica, MA · USA  
Phone +1 (978) 439-9899  
Fax +1 (978) 663-9177  
info@brukeroptics.com

### Bruker Shanghai Ltd.

Shanghai · China  
Phone +86 21 51720-890  
Fax +86 21 51720-899  
info.bopt.cn@bruker.com